

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RAC1 Human, Recombinant, E.coli

产品货号: TP03653

第三版

别名: Ras-related C3 botulinum toxin substrate 1 isoform Rac1, MIG5, p21-Rac1, TC-25

描述: Rac1, belonging to the Rho small GTPase family, regulates the actin cytoskeleton but also other cellular processes. RAC1 have been shown to be involved in the regulation of cell-cell adhesion. Recombinant RAC1 was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In 20 mM Tris-HCl buffer (pH 7.5) containing 1 mM DTT 2 mM EDTA

分子量: 21.4 kDa (192 aa), confirmed by MALDI-TOF

序列:

MQAIKCVVVG DGAVGKTC LLISYTTNAFPGEYIPTVFDNYSANVMVDGKPVNLGLWDTAGQEDYDRLRPLSYPQTDVFL
ICFSLVSPASFENVRAKWYPEVRHHCPNTPHILVGTKLDRDDKDTIEKLKEKCLTPITYPQGLAMAKEIGAVKYLECSALTQ
RGLKTVFDEAIRAVLCPPPVKRKRKCLLL

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: < 1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.