

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

RALY, 1-306aa, Human, His tag, E.coli

产品货号: TP03666

第三版

别名:RNA-binding protein Raly, RNA binding protein Raly, HNRPCL2, P542

描述:RALY is a member of the heterogeneous nuclear ribonucleoprotein (hnRNP) gene family. This protein may play a role in pre-mRNA splicing and in embryonic development. Alternate splicing results in multiple transcript variants. Recombinant human RALY protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

分子量:34.9 kDa (329aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMCLKLQASNVTKNDPKSINSRVFIGNLNTALVKKSDVETIFSKYGRVAGCSVHKG
YAFVQYSNERHARAAVLGENGRVLAGQTLINMAGEPKPDRPKGLKRAASAIYSGYIFDYDYRDFYDRLFDYRGRSLP
VPVPRAVPVKRPRVTVPLVRRVKTNPVKLFARSTAVTTSSAKIKLKSELQAIKTELQIKSNIDALLSRLEQIAAEQKANP
DGKKKGDGGGAGGGGGGGGGGGGGSSRPPAPQENTTSEAGLPQGEARTRDDGDEEGLLTHSEEELEHSQ
DTDADDGALQ

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.