



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Rantes, 24-91 aa, Human, Recombinant, E.coli

产品货号: TP03673

第三版

别名: Small inducible cytokine A5, CCL5, D17S136E, SCYA5, SISd, TCP228.

描述: Rantes, also known as CCL5, is the protein classified as a chemotactic cytokine or chemokine. This protein is a chemotactic for blood monocytes, memory T helper cells, eosinophils and basophils, and plays an active role in recruiting leukocytes into inflammatory sites. Rantes also induces the proliferation and activation of certain natural-killer (NK) cells to form CHAK (CC-Chemokine-activated killer) cells.

配方: Liquid. In 10 mM Sodium Citrate (pH 3.5) containing 10% Glycerol.

分子量: 7.9 kDa (69aa)

序列: MSPYSSDTTPCCFAYIARPLPRAHIKEYFYTSGKCSNPAVVVFVTRKNRQVCANPEKKWVREYINSLEMS

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.