

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

RASL12, 1-266aa Human, His tag, E.coli

产品货号: TP03684

第三版

**别名:**Ras-like protein family member 12, RIS, Ras-like protein Ris, Ras-like family 12

**描述:**RASL12 belongs to the small GTPase superfamily, Ras family. Recombinant human RASL12 protein, fused to His-tag at N-terminus, was expressed in E.coli.

**配方:**Liquid, In 20mM Tris-HCl (pH8.0) containing 10% glycerol

**分子量:**32.1 kDa (289aa)

**序列:**

MGSSHHHHHHSSGLVPRGSHMGSMSSVFGKPRAGSGPQSAPLEVNLAAILGRRGAGKSALTVKFLTKRFISEYDPNLEDTYS  
SEETVDHQPVHLRVMDTADLDTPRNCERYLNWAHAFLVVYSVDSRQSFSSSYLELLALHAKETQRSIPALLLGNKLDL  
AQYRQVTKAEGVALAGRFGLFFFEVSACLD FEHVQHFHEAVREARRELEKSPLTRPLFISEERALPHQAPL TARHGLASC  
TFNTLSTINLKEMPTVAQAKLVTVKSSRAQSKRKAPTLLKGFKIF

**纯度:**> 95% by HPLC

**浓度:**0.5 mg/ml (determined by Bradford assay)

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.