

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RASSF3, 1-238aa, Human, His tag, E.coli

产品货号: TP03686

第三版

**别名:**Ras association domain-containing protein 3, MGC119194, MGC119195, MGC119197, Ras association (RalGDS/AF 6) domain family 3 antibody

**描述:**Ras association domain-containing protein 3, also known as RASSF3, is a member of RASSF family. Members of the RASSF family contain a highly conserved Ras association domain (Ral GDS/AF-6 or RA) and function as Ras effectors/tumor suppressors. RASSF3 is a ubiquitously expressed protein found in normal and cancerous tissues. It contains an N-terminal RA domain and a coiled-coil SARAH domain. Three isoforms, namely RASSF3A, RASSF3B and RASSF3C, may exist for RASSF3 due to alternative splicing. RASSF3B and RASSF3C are shorter than RASSF3A and do not contain the RA and SARAH domain. Recombinant human RASSF3 protein, fused to His-tag at N-terminus, was expressed in E.coli.

**配方:**Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.4M Urea, 10% glycerol

**分子量:**30 kDa (261aa)

**序列:**

MGSSHHHHHSSGLVPRGSHMGSMSSGYSSLEEDAEDFFFTARTSFFRRAPQGGKPRSGQQDVEKEKETHSYLSKEEIKEKV  
HKYNLAVTDKMKMTLNSNGIYTGFIKVMELCKPPQTSPNSGKLSPPSSNGCMNTLHISSTNTVGEVIEALLKKFLVTESPAK  
FALYKRCHREDQVYACKLSDREHPLYLRLVAGPRTDTLSFVLRHEHEIGEWAFSLPELQNFRLILDKEEDEQLQNLKRRRYTA  
YRQKLEELREVWVKPD

**纯度:**> 95% by HPLC

**浓度:**1 mg/ml (determined by Bradford assay)

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:**+4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.