

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RBBP4, 1-425aa Human, His tag, E.coli

产品货号: TP03687

第三版

别名:Retinoblastoma binding protein 4, NURF55, RBAP48.

描述:Retinoblastoma binding protein 4, also known as RBBP4, is a ubiquitously expressed nuclear protein which belongs to a highly conserved subfamily of WD-repeat proteins. It is present in protein complexes involved in histone acetylation and chromatin assembly. RBBP4 is also part of co-repressor complexes, which is an integral component of transcriptional silencing.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% glycerol, 1mM EDTA

分子量:50.2 kDa (449aa)

序列:

MGSSHHHHHSSGLVPRGSHMGSHMADKEAAFDDAVEERVINEEYKIWKKNTPFLYDLVMTHALEWPSLTAQWLDPV
TRPEGKDFSIHRLVLGHTSDEQNHLVIASVQLPNDDAQFDASHYDSEKGEFGGFGSVSGKIEIEIKINHEGEVNRARYMP
QNPCIATKTPSSDVLVFDYTKHPSKPDPSGECNPDRLRGRHQKEGYGLSWNPNLSGHLLSASDDHTICLWDISAVPKEGK
VVDAKTIFTGHTAVVEDVSWHLLHESLFGSVADDQKLMIWDRSNNTSKPSHSVDAHTAEVNCLSFNPYSEFILATGSAD
KTVALWDLRNLKLLHSFESHKDEIFQVQWSPHNETILASSGTDRLNVWDLISKIGEEQSPEDAEDGPPELLFIHGGHTAK
ISDFSWNPNEPWVICSVSEDNIMQVWQMAENIYNDEDPEGSVDPEGQGS

纯度:> 95% by HPLC

浓度:0.5 mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.