



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RBM18, 1-190aa, Human, His tag, E.coli

产品货号: TP03692

第三版

别名: Probable RNA-binding protein 18, MGC2734, RNA-binding motif protein 18, RNA binding motif protein 18

描述: RBM18 contains 1 RRM (RNA recognition motif) domain. Recombinant human RBM18 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography.

配方: Liquid. Phosphate Buffered Saline (pH7.4) containing 10% glycerol

分子量: 24.0 kDa(213aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMEAETKTLPLENASILSEGLQEGHRLWIGNLDPKITEYHLLKLLQKFGKVKQFDLFLF  
HKSGALEGQPRGYCFVNFETKQEAQAIQCLNGKLALS KLVVRWAHAQVKRYDHNKNDKILPISLEPSSSTEPTQSNLSV  
TAKIKAIEAKLKMMMAENPDAEYPAAPVYSYFKPPDKKRTTPYSRTAWKSRR

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.