

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

RNMTL1, 1-420aa, Human, His tag, E.coli

产品货号: TP03760

第三版

别名:RNA methyltransferase-like protein 1, HC90

描述:RNMTL1 also known as RNA methyltransferase-like protein 1 is belong to the RNA methyltransferase TrmH family. This protein is probable RNA methyltransferase and expressed at same level in normal liver and hepatocarcinoma. Recombinant human RNMTL1, fused to His-tag at N-terminus, was expressed in E.coli.

配方:Liquid, In 20mM Tris-HCl (pH8.0) containing 10% glycerol

分子量:49.4kDa (443aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMALVPRARFVVRPLLQVVQAWDL DARRWVRALRRSPVKVVFPSGEVVEQKRAP
GKQPRKAPSEASAQEQREKQPLEESASRAPSTWEESGLRYDKAYPGDRRLSSVMTIVKSRPFREKQGKILLEGRRLISDALK
GAVPKMFFFSRLEYLKELPVDKLGVS LKVKFEDIKDWSDLVTPQGIMGIFAKPDHVKMTYPKTQLQHSLPLLLICDNL RD
PGNLGTILRSAAGAGCSKVL LTKGCVD AWE PKVLRAGMGAHFRMPIINNLEWETV PNYLPPDTRVYVADNCGLYAQAE
MSNKASDHGWVCDQRVMKFHKYEEEEEDVETGASQDWLPHVEVQSYDS DWTEAPAAVVIGGETYGV SLES LQLAESTG
GKROLLIPVVPVGD SLNSAMAASILLFEGKRQLRGRAEDLSRDRSYH

纯度:> 95% by HPLC

浓度:1mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.