

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RPL11, 1-178aa, Human, His tag, E.coli

产品货号: TP03767

第三版

别名: Ribosomal protein L11, DBA7, GIG34, L11

描述: Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. RPL11 is a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L5P family of ribosomal proteins. It is located in the cytoplasm. The protein probably associates with the 5S rRNA. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Recombinant human RPL11 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% glycerol, 1mM DTT, 2mM EDTA, 250mM Imidazole

分子量: 22.6kDa (201aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHSSGLVPRGSHMGSMADQGEKENPMRELRIRKLCNICVGESGDRLTRAAKVLEQLTGQTPVFSKARY
TVRSFGIRRNEKIAVHCTVRGAKAEIILEKGLKVREYELRKNFSDTGNFGFGIQEHIDLGIKYDPSIGIYGLDFYVVLGRPG
FSIADKKRRTGCIGAKHRISKEEAMRWFQQKYDGIILPGK

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by BRADFORD assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.