

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RPL12, 1-165aa, Human, His tag, E.coli

产品货号: TP03768

第三版

别名: Ribosomal protein L12, L12

**描述:** Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. RPL12 is a ribosomal protein that is a component of the 60S subunit. This protein belongs to the L11P family of ribosomal proteins. It is located in the cytoplasm. RPL12 binds directly to the 26S rRNA. Recombinant human RPL12 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

**配方:** Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1mM DTT

**分子量:** 20.2 kDa (188aa) confirmed by MALDI-TOF

**序列:**

MGSSHHHHHSSGLVPRGSHMGSMPPKFDPNKIKVVYLRCTGGEVVGATTSALAPKIGPLGLSPKKVGGDDIAKATGDWKGL  
RITVKLTIQNRQAQIEVVPASALIIKALKEPPRDRKKQKNIKHSGNITFDEIVNIARQMRHRSLARELSGTIKEILGTAQSVG  
CNVDGRHPHDIIDDINSGAVECPAS

**纯度:** > 95% by HPLC

**浓度:** 0.5 mg/ml (determined by Bradford assay)

**内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.