

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

RPL26L1, 1-145aa, Human, His tag, E.coli

产品货号: TP03772

第三版

别名:60S ribosomal protein L26-like 1, PL26P1.

描述:60S ribosomal protein L26-like 1, also known as RPL26L1, belongs to the ribosomal protein family and may play a role in protein synthesis. It is located in the cytoplasm. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 40% glycerol, 2mM DTT

分子量:19.6 kDa (168aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMKFNPFVTSRDRSKNRKRHFNAPSHVRRKIMSSPLSKELRQKYNVRSMPIRKDDEVQ
VVRGHYKGGQIGKVVQVYRKKYVIYIERVQREKANGTTVHVGIHPSKVVITRLKLDKDRKKILERKAKSRQVGKEKGY
KEELIEKMQE

纯度:> 95% by HPLC

浓度:0.25 mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.