

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

RPS20, 1-142aa, Human, His tag, E.coli

产品货号: TP03793

第三版

别名:40S ribosomal protein S20 isoform 1, S20

描述:Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S10P family of ribosomal proteins. It is located in the cytoplasm. This gene is co-transcribed with the small nucleolar RNA gene U54, which is located in its second intron. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Two transcript variants encoding different isoforms have been identified for this gene. Recombinant human RPS20 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% glycerol, 2mM DTT, 1mM EDTA

分子量:18.4kDa (165aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHSSGLVPRGSHMGSMFAFKDTGKTPVEPEVAIHRIRITLTSRNVKSLEKVCADLIRGAKEKNLKVKGPPVRM
PTKTLRITTRKTPCGEGSKTWDRFQMRIHKRLIDLHSPSEIVKQITSISIEPGVELIESTDAEPMDEGQQYTLRSVFESPGTCP
F

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.