

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RPS4X, 1-263aa, Human, His tag, E.coli

产品货号: TP03797

第三版

别名: 40S ribosomal protein S4, X isoform X isoform, Ribosomal protein S4, X-linked, CCG2; DXS306, RPS4, S4, SCAR, SCR10

**描述:** Cytoplasmic ribosomes, organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. RPS4X is ribosomal protein S4, a component of the 40S subunit. Ribosomal protein S4 is the only ribosomal protein known to be encoded by more than one gene, namely this gene and ribosomal protein S4, Y-linked (RPS4Y). The 2 isoforms encoded by these genes are not identical, but are functionally equivalent. Ribosomal protein S4 belongs to the S4E family of ribosomal proteins. This gene is not subject to X-inactivation. It has been suggested that haploinsufficiency of the ribosomal protein S4 genes plays a role in Turner syndrome; however, this hypothesis is controversial. Recombinant human RPS4X protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

**配方:** Liquid. In 20mM Tris-HCl(pH8.0) containing 40% glycerol, 0.15M NaCl, 1mM DTT

**分子量:** 32kDa (286aa), confirmed by MALDI-TOF

**序列:**

MGSSHHHHHSSGLVPRGSHMGSMARGPKKHLKRVAAPKHWMLDKLTGVFAPRPSTGPHKLRECLPLIIFLRNRLKYAL  
TGDEVKKICMQRFIKIDGKVRTDITYPAGFMDVISIDKTGENFRLIYDTKGRFAVHRITPEEAKYKLCVKRIFVGTGKIPHL  
VTHDARTIRYPDPLIKVNDTIQIDLETGKITDFIKFDTGNLCMVTGGANLGRIGVITNRERHPGSFDVVHVKDANGNSFAT  
RLSNIFVIGKGNKPWISLPRGKGIRLTIAEERDKRLAAKQSSG

**纯度:** > 95% by HPLC

**浓度:** 0.5 mg/ml (determined by Bradford assay)

**内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.