

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

S100A1, 1-94aa, Mouse, His tag, E.coli

产品货号: TP03817

第三版

别名:S100 calcium binding protein A1, S100, S100a

描述:S100 calcium binding protein A1, also known as S100A, is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100A1 is expressed by striated muscle, endothelium and chondrocytes.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 30% glycerol, 0.1M NaCl

分子量:12.6 kDa (114aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSELESAMETLINVFHAHSGKEGDKYKLSKKELKDLLQTELSGFLDVQKDADAVDKV
MKELDENGDFKEYVVLVAALTVACNNFFWETS

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储:+4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.