



本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

SERF2, 1-59aa, Human, His tag, E.coli

产品货号: TP03894

第三版

别名:Small EDRK-rich factor 2, 4F5REL FAM2C, H4F5REL, HsT17089

描述:SERF2, also known as small EDRK-rich factor2, belongs to the SERF family. Recombinant human SERF2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 30% glycerol

分子量:9.3 kDa (82aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMTRGNQRELARQKNMKKQSDSVK GKRRDDGLSAAARKQRDSEIMQQKQKKANE
KKEEPK

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by BCA assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.