

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

SH3BGRL2, 1-107aa, Human, His tag, E.coli

产品货号: TP03924

第三版

别名:SH3 domain-binding glutamic acid-rich-like protein 2, FASH3, Fovea-associated SH3 domain-binding protein, SH3 domain-binding glutamic acid-rich-like protein 2, SH3 domain binding glutamic acid rich like protein 2.

描述:SH3 domain-binding glutamic acid-rich-like protein 2, also known as SH3BGRL2, is a member of the human SH3BGR family. SH3BGRL2 gene located to chromosome 6q13-15 and its messenger RNA has a large 3' untranslated region containing several AUUUA repeats. SH3BGRL2 codes for a protein of 107 amino acids, which, like SH3BGRL and SH3BGRL3 proteins, is highly homologous to the N-terminal region of the SH3BGR protein and appears to be related to Glutaredoxins and to PKC-interacting cousin of thioredoxin homology domain. Recombinant human SH3BGRL2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1mM DTT

分子量:14.7 kDa (130aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHSSGLVPRGSHMGSMVIRVFIASSSGFVAIKKKQQDVVRFLEANKIEFEEVDITMSEEQRQWMYKNVPPE
KKPTQGNPLPPQIFNGDRYCGDYDSFFESKESNTVFSFLGLKPRLASKAEP

纯度:> 95% by HPLC

浓度:0.5 mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.