

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

SLA2, 1-261aa, Human, His tag, E.coli

产品货号: TP03955

第三版

别名:Src-like-adaptor 2 isoform a, Src-like-adaptor 2, C20orf156, MARS, SLAP-2, SLAP2

描述:SLA2 is a member of the SLAP family of adapter proteins. The encoded protein may play an important receptor-proximal role in downregulating T and B cell-mediated responses and inhibits antigen receptor-induced calcium mobilization. This protein interacts with Cas-Br-M (murine) ecotropic retroviral transforming sequence c. Two transcript variants encoding distinct isoforms have been identified for this gene. Recombinant human SLA2 protein, fused to His-tag at N-terminus, was expressed in E.coli .

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing , 10% glycerol, 0.4M Urea

分子量:31.0kDa (284aa)

序列:

```
MGSSHHHHHHSSGLVPRGSHMGSMGSLPSRRKSLPSPSLSSSVQGGPVTMEAERSKATAVALGSFPAGGPAELSLRLGEP  
LTIVSEGDWWTVLSEVSGREYNIPSVHVAKVSHGWLYEGLSREKAEELLLLPGNPGGAFLIRESQTRRGSYSLSVRLSRPAS  
WDRIRHYRIHCLDNGWLYISPRLTFPSLQALVDHYSELADDICLLKEPCVLQRAGPLPGKDIPLPVTVQRTPLNWKELDSS  
LLFSEAAATGEESLLSEGLRESLSFYISLNDEAVSLDDA
```

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.