

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

SNRPD2, 1-118aa, Human, His tag, E.coli

产品货号: TP03988

第三版

别名: Small nuclear ribonucleoprotein polypeptide D2, Sm-D2, SMD2, SNRPD1

描述: SNRPD2 is a small nuclear ribonucleoprotein (snRNPs) that comprises the spliceosome in eukaryotes. It is required for pre-mRNA splicing and small nuclear ribonucleoprotein biogenesis. Alternative splicing occurs at this locus and two transcript variants encoding the same protein have been identified. Recombinant human SNRPD2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography.

配方: Liquid. 20mM Tris-HCl buffer (pH8.0) containing 50% glycerol 0.2M NaCl, 1mM DTT.

分子量: 15.6 kDa(138aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHSSGLVPRGSHMSLLNPKKSEMTPEELQKREEEFNTPGLSVLTQSVKNNTQVLINCRNNKLLGRVKAF
DRHCNMVLENVKEMWTEVPKSGKGGKKKSKPVNKDRYISKMFLRGDSVIVVLRNPLIAGK

纯度: > 95% by HPLC

浓度: 0.25 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.