

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

SNURF, 1-71aa, Human, His tag, E.coli

产品货号: TP03996

第三版

别名:SNRPN upstream reading frame protein, Small nuclear ring finger protein

描述:SNURF is a highly basic protein localized to the nucleus. The evolutionarily constrained open reading frame is found on a bicistronic transcript which has a downstream ORF encoding the small nuclear ribonucleoprotein polypeptide N. The upstream coding region utilizes the first three exons of the transcript, a region that has been identified as an imprinting center. Multiple transcription initiation sites have been identified and extensive alternative splicing occurs in the 5' untranslated region but the full-length nature of these transcripts has not been determined. An alternate exon has been identified that substitutes for exon 4 and leads to a truncated, monocistronic transcript. Alternative splicing or deletion caused by a translocation event in the 5' UTR or coding region of this gene leads to Angelman syndrome or Prader-Willi syndrome due to parental imprint switch failure. The function of this protein is not yet known. Recombinant human SNURF protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 40% glycerol, 1mM DTT.

分子量:10.8 kDa (94aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHSSGLVPRGSHMGSMERARDRLHLRRTTEQHVPEVEVQVKRRRTASLSNQECQLYPRRSQQQVPVVD
FQAELRQAFLAETPRGG

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.