

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

SPINK7, 20-85aa, Human, His tag, E.coli

产品货号: TP04016

第三版

别名: Serine peptidase inhibitor Kazal type 7, ECG2, ECRG2

描述: Serine peptidase inhibitor Kazal type 7, also known as SPINK7, contains 1 Kazal-like domain and is probably serine protease inhibitor. Recombinant human SPINK7 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: In 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 40% glycerol, 2mM DTT, 0.1mM PMSF, 50mM imidazole

分子量: 9.6 kDa (89aa), confirmed by MALDI-TOF

序列:

纯度: > 95% by HPLC

浓度: 0.25 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.