

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

SRS (Seryl-tRNA synthetase)

产品货号: TP04034

第三版

别名:

**描述:**Seryl-tRNA synthetase catalyzes the ligation of serine with seryl-tRNA using ATP. This gene belongs to the class II amino-acyl tRNA family. The encoded enzyme catalyzes the transfer of L-serine to tRNA (Ser) and is related to bacterial and yeast counterparts. Multiple alternatively spliced transcript variants have been described but the biological validity of all variants is unknown.

**配方:**1X PBS (pH7.2) 200mM NaCl 20% Glycerol

**分子量:**58kDa

序列:

**纯度:**> 95% by HPLC

**浓度:**0.50mg/ml

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.