

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

SSSCA1, 1-199aa, Human, His tag, E.coli

产品货号: TP04040

第三版

别名:Sjogren syndrome/scleroderma autoantigen 1, p27

描述:SSSCA1 is recognized by a subset of anti-centromere antibodies from patients with scleroderma and/or Sjogren's syndrome. Subcellular localization has not yet been established. Recombinant human SSSCA1 protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方:Liquid. 20mM Tris-HCl buffer (pH8.0) containing 10% glycerol 0.4M Urea

分子量:23.9kDa (222aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMALNGAEVDDFSWEPTEAETKVLQARRERQDRISRLMGDYLLRGYRMLGETCA  
DCGTILLQDKQRKIYCVACQELSDVDKDNPALNAQAALSQAREHQLASASELPLGSRPAPQPPVPRPEHCEGAAAGLKA  
AQGPPAPAVPPNTDVMACTQTALLQKLTWASAELGSSTSLETSIQLCGLIRACAEALRSLQQLQH

纯度:> 95% by HPLC

浓度:0.5mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.