

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

STOML1, 79-398 aa, Human, His tag, E.coli

产品货号: TP04070

第三版

别名: Stomatin-like protein 1, huNC-24, SLP-1, STORP

描述: STOML1 belongs to the band 7/mec-2 family and contains 1 SCP2 domain. This protein is ubiquitously expressed at low levels. Recombinant human STOML1 protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方: Liquid. In 20mM Tris-HCl (pH8.0) containing 10% glycerol

分子量: 37.0kDa (343aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSLKIVPTYERMIVFRLGRIRTPQGPGMVLLLLPFIDSFQRVDLRTRAFNVPPCKLASKDG  
AVLSVGADVQFRIWDPVLSVMTVKDLNTATRMTAQNAMTKALLKRPLREIQMEKLIKISDQLLEINDVTRAWGLEVDV  
ELAVEAVLQPPQDSPAGPNLDSTLQQLALHFLGGSMNSMAGGAPSPGPADTVEMVSEVEPPAPQVVGARSSPKQPLAEGLL  
TALQPFLSEALVSQVGACYQFNVLPSGTQSAYFLDLTTGRGRVGHGVPDGPVVEVMAEADLRALLCRELRPLGAYMS  
GRLKVKGDLAMAMKLEAVLRALK

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.