

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

STX11, 1-287aa, Human, His tag, E.coli

产品货号: TP04074

第三版

别名:Syntaxin-11, FHL4, HLH4, HPLH4

描述:STX11, also known as syntaxin-11, is a member of the t-SNARE family. It acts to regulate protein transport between late endosomes and the trans-Golgi network. Recombinant human STX11 protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方:Liquid. 20mM Tris-HCl buffer (pH8.0) containing 10% glycerol, 0.4M Urea

分子量:35.8 kDa (311aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMKDRLAELLDLSKQYDQQFPDGDDEFDSPHEDIVFETDHILESLYRDIRDIQDEN
QLLVADV KRLGKQNRFLTSMRRLSSIKRDTNSIAKAIKARGEVIHCKLRAMKELSEAAEAQHGPSAVARISRAQYNALT
LTFQRAMHDYNQAEMKQRDNCKIRIQRQLEIMGKEVSGDQIEDMFEQGWVDFSENLLADVKGARAALNEIESRHREL
LRLESRI RDVHELFLQMAVLVEKQADTLNVIELNVQKTVDYTGQAKAQVRKAVQYEEKNPCRTLCCFCCPCLK

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.