

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

SULT2B1, 1-365aa, Human, E.coli

产品货号: TP04093

第三版

别名: HSST2, ST2B1, Sulfotransferase family cytosolic 2B member 1 isoform b Alcohol sulfotransferase, EC 2.8.2.2, Hydroxysteroid sulfotransferase 2, Sulfotransferase 2B1, Sulfotransferase family cytosolic 2B member 1, SULT2B1a, SULT2B1b.

描述: Sulfotransferase family cytosolic 2B member 1 isoform b, also known as SULT2B1, is a member of Sulfotransferase family. Sulfotransferase enzymes catalyze the sulfate conjugation of many hormones, neurotransmitters, drugs, and xenobiotic compounds. These cytosolic enzymes are different in their tissue distributions and substrate specificities. SULT2B1b is localized in the cytosol and nuclei of human cells. SULT2B1b is selective for the sulfation of 3beta-hydroxysteroids such as dehydroepiandrosterone and pregnenolone, and may also have a role in cholesterol sulfation in human skin. Recombinant human SULT2B1, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In 20mM Tris-HCl buffer (pH7.5) containing 10% glycerol

分子量: 41.3 kDa (365aa), confirmed by MALDI-TOF.

序列:

MDGPAEPQIPGLWDTYEDDISEISQKLPGEYFRYKGVPPFVGLYSLESISLAENTQDVRDDDFIITYPKSGTTWMIIEICLILK
EGDPSWIRSVPIWERAPWCETIVGAFSLPDQYSPRLMSSHLPQIFTKAFFSSKAKVIYMGRRNPRDVVVSLYHYSKIAGQLKD
PGTPDQFLRDFLKGVEVQFGSWFDHIKGWLRMKGKDNFLFITYEELQQDLQGSVERICGFLGRPLGKEALGSVVAHSTFSA
MKANTMSNYTLLPSSLDDHRRGAFLRKGVCGDWKNHFTVAQSEAFDRAYRKQMRGMPTFPWDEDPEEDGSPDPEPSPE
PEPKPSLEPNTSLEREPNNSPSPSPGQASETPHPRPS

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.