

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

SYCE3, 1-88aa, Human, His tag, E.coli

产品货号: TP04105

第三版

别名:synaptonemal complex central element protein 3, testis highly expressed protein 2, C22orf41, chromosome 22 open reading frame 41, THEG-2, THEG2, SYCE-3

描述:SYCE3(synaptonemal complex central element protein 3) is major component of the transverse central element of synaptonemal complexes (SCS). The synaptonemal complex is a protein structure that forms between homologous chromosomes (two pairs of sister chromatids) during meiosis and that is thought to mediate chromosome pairing, synapsis, and recombination (crossing-over). This protein required for chromosome loading of the central element-specific SCS proteins, and for initiating synapsis between homologous chromosomes. Recombinant human SYCE3 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 20% glycerol, 2mM DTT

分子量:12.8 kDa (108aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHSSGLVPRGSHMDDADPEERNYDNMLKMLSDLNKDLEKLLLEEMEKISVQATWMAYDMVVMRTNPTLA
ESMRRLEDAFVNCKEEMKKNWQELLHETKQRL

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.