

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

-Synuclein 1-60, Human, Recombinant, E.coli

产品货号: TP04112

第三版

别名: SNCA, NACP, PARK1, alpha-Synuclein, Non-A beta component of AD amyloid, Non-A4 component of amyloid precursor, Alpha synuclein, Alpha-synuclein isoform NACP140, alphaSYN, MGC105443, MGC110988, MGC127560, MGC64356

描述: A deletion mutant of alpha-synuclein (amino acids 1-60), which contains the N-terminal amphipathic domain. alpha-Syn 1-60 was overexpressed in E. coli and the recombinant protein was purified to apparent homogeneity by using conventional column chromatography techniques.

配方: Liquid. In 20 mM Tris-HCl buffer (pH 7.5) containing 0.1 M NaCl

分子量: 6.1 kDa (60 aa), confirmed by MALDI-TOF.

序列:

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by BCA assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.