

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

-Synuclein -NAC, Human, Recombinant, E.coli

产品货号: TP04121

第三版

别名:SNCA, NACP, PARK1, alpha-Synuclein Delta-NAC, Non-A beta component of AD amyloid, Non-A4 component of amyloid precursor, NACP, Alpha-synuclein, Alpha synuclein, Alpha-synuclein, isoform NACP140, alphaSYN, MGC105443

描述:A deletion mutant of alpha-synuclein that lacks the NAC region(amino acid 61-95). alpha-synuclein (deltaNAC) was overexpressed in E. coli and purified to apparent homogeneity by taking advantage of the thermosolubility of the protein and by using conventional column chromatography techniques. 6 amino acids were added as a linker.

配方:Liquid. In 20 mM Tris-HCl buffer (pH 7.5) containing 0.1 M NaCl

分子量:11.9 kDa (111 aa), confirmed by MALDI-TOF. (Molecular weight on SDS-PAGE will appear higher)

序列:

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by BCA assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.