

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

TACO1, 61-297aa, Human, His tag, E.coli

产品货号: TP04132

第三版

别名: Translational activator of cytochrome c oxidase 1, CCDC44

描述: TACO1 is a mitochondrial protein that function as a translational activator of mitochondrially-encoded cytochrome c oxidase 1. Mutations in this gene are associated with Leigh syndrome. Recombinant human TACO1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1mM DTT

分子量: 28.8 kDa (260aa) confirmed by MALDI-TOF

序列:

纯度: > 95% by HPLC

浓度: 0.5mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.