

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

TCEAL8, 1-117aa , Human, His tag, E.coli

产品货号: TP04166

第三版

别名:transcription elongation factor A (SII)-like 8,

描述:TCEAL8 (Transcription elongation factor A protein-like 8), also known as TCEA-like protein 8, belongs to the TFS-II family and TFA subfamily. This protein is involved in transcriptional regulation and localized in nucleus. Recombinant human TCEAL8 protein, fused to His-tag at C-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. In 20 mM Tris-HCl buffer, pH8.0, 20% glycerol, 1mM DTT, 200mM NaCl

分子量:14.7kDa (125aa), confirmed by MALDI-TOF (Molecular weight on SDS-PAGE will appear higher)

序列:

MQKSCEENEGKPKQNMPKAEEDRPLEDVPQEAEGNPQPSEEGVSQEAEGNPRGGPNQPGQGFKEDTPVRHLDPEEMIRG  
VDELERLRREEIRRVRNKFMHMMHWKQRHSRSPYPVCFRPLEHHHHHHH

纯度:> 95% by HPLC

浓度:0.5mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.