

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

TPRKB, 1-175 6aa, Human, His tag, E.coli

产品货号: TP04290

第三版

别名:TP53RK binding protein, CGI-121

描述:TPRKB belongs to the CGI121/TPRKB family. It is widely expressed and localized to nucleus and cytoplasm. This protein Interacts with TP53RK/PRPK. Recombinant human TPRKB protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. 20mM Tris-HCl buffer (pH8.0) containing 10% glycerol 0.15M NaCl

分子量:22.2kDa(199aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMQLTHQLDLFPECRVTLTLLFKDVKNAGDLRRKAMEGTIDGSLINPTVIVDPFQIL
VAANKAVHLYKLGKMKTRTLSTEIIFNLSPNNNISEALKKFGISANDTSILIVYIEEGEKQINQEY LISQVEGHQVSLKLNLP
EIMNITEVKKIYKLSSQEESIGTLLDAIICRMSTKDVL

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.