

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

TRS (Threonyl-tRNA synthetase)

产品货号: TP04315

第三版

别名:

描述:Aminoacyl-tRNA synthetases catalyze the aminoacylation of tRNA by their cognate amino acid. Because of their central role in linking amino acids with nucleotide triplets contained in tRNAs, aminoacyl-tRNA synthetases are thought to be among the first proteins that appeared in evolution. Threonyl-tRNA synthetase belongs to the class-II aminoacyl-tRNA synthetase.

配方:1X PBS (pH7.2) 200mM NaCl 20% Glycerol

分子量:83kDa

序列:

纯度:> 95% by HPLC

浓度:0.50mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.