

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

TSR2, 1-190aa, Human, His tag, E.coli

产品货号: TP04327

第三版

别名: Regulation of nuclear pre-mRNA domain-containing protein 1A, DT1P1A10, RP1-112K5.2, WGG1

描述: TSR2 belongs to the TSR2 family. This protein may be involved in 20S pre-rRNA processing. Recombinant human TSR2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In PBS buffer (pH 7.4) containing 10% glycerol

分子量: 23.2kDa (213aa) confirmed by MALDI-TOF (Molecular size on SDS-PAGE will appear higher)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMAGAAEDARALFRAGVCAALEAWPALQIAVENGFGGVHSQEKAKWLGGAVEDY  
FMRNADLELDEVEDFLGELLTNEFDTVVEDGSLPQVSQQQLQTMFHHFQRGDGAALREMACITQRKCKVTATALKTARE  
TDEDEDDVDSVEEMVATNDGAATDGVCPEPSPDAQTIKEEDIVEDGWTIVRRK

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.