



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

TUSC2, 1-110aa, Human, His tag, E.coli

产品货号: TP04352

第三版

别名:Tumor suppressor candidate 2, Tumor suppressor candidate 2, C3orf11, FUS1, PAP, PDAP2

描述:TUSC2 belongs to the TUSC2 family. This protein may function as a tumor suppressor, inhibiting colony formation, causing G1 arrest and ultimately inducing apoptosis in homozygous 3p21.3 120-kb region-deficient cells. Recombinant human TUSC2 protein, fused to His-tag at N-terminus, was expressed in E.coli .

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.4M Urea

分子量:14.5 kDa (133aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMGASGSKARGLWPFASAAGGGGSEAAGAEQALVRPRGRAVPPFVFTRRGSMFYDE
DGDLAHEFYEETIVTKNGQKRAKLRRVHKNLIPQGIVKLDHPRIHVDFPVILYEV

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.