

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Ubc7(Ubiquitin Conjugating enzyme 7) Human Recombinant, expressed in E.coli

产品货号: TP04372

第三版

别名:UBE2G1, UBE2G, Ubc7 (Ubiquitin Conjugating enzyme 7), EC 6.3.2.19, Ubiquitin-protein ligase G1, Ubiquitin carrier protein G1, E217K, UBC7, Ubiquitin-conjugating enzyme E2 G1, Ubiquitin-conjugating enzyme E2L 3 UBC 7

描述:Human Ubiquitin-conjugating enzyme 7(UbcH7) is a ubiquitin-conjugating enzyme (E2) mediating c-fos degradation, transcription factor NF-kappa B maturation, and human papilloma virus-mediated p53 and Myc protein degradation, in vitro. The ubiquitin-conjugating enzymes(E2s) are essential components of the post-translational protein ubiquitination pathway, mediating the transfer of activated ubiquitin to substrate proteins. The human UBE2L1-UBE2L4 gene could potentially encode different isoforms of the UbcH7. UBE2L3 gene, located at chromosome 22q11.2, is the only identical family member with introns and encodes a polypeptide sequence identical to that of UbcH7. UbcH7(154amino acid) was over-expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. In 50 mM HEPES (pH7.5) 150 mM NaCl, 1mM DTT, 10% glycerol

分子量:17.9 kDa (154 aa), confirmed by MALDI-TOF.

序列:

MAASRRLMKELEEIRKCGMKNFRNIQVDEANLLTWQGLIVPDNPPYDKGAFRIEINFPAEYPFKPPKITFKTKIYHPNIDEK
GQVCLPVISAENWKPATKTDQVIQSLIALVNDPQPEHPLRADLAEYISKDRKKFCKNAEEFTKKYGEKRPVD

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.