

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Ubc9 (Human ubiquitin conjugating enzyme 9) Human, Recombinant, expressed in E.coli

产品货号: TP04373

第三版

别名:UBE2I, UBC9, UBCE9, Ubc9, Human ubiquitin conjugating enzyme 9, EC 6.3.2, SUMO-protein ligase, Ubiquitin-conjugating enzyme E2 I, Ubiquitin-protein ligase I, Ubiquitin carrier protein I

描述:Human Ubc9 is homologous to ubiquitin-conjugating enzymes(E2s). However, instead of conjugating ubiquitin, it conjugates a ubiquitin homologue, small ubiquitin-like modifier 1(SUMO-1). And hUbc9 retains striking structural and functional conservation with yeast Ubc9. The ubiquitin-dependent protein degradation system has been recognized as a complete enzymatic pathway that is responsible for the selective degradation of abnormal and short-lived proteins. The conjugation of ubiquitin requires the activities of ubiquitin-activating(E1) and -conjugating(E2) enzymes.

配方:Liquid. In 50 mM HEPES (pH7.5) 150mM NaCl, 1mM DTT, 10%glycerol.

分子量:18 kDa (158 aa), confirmed by MALDI-TOF.

序列:

MSGIALSRLAQERKAWRKDHPFGFVAVPTKNPDGTMNLMNWECAIPGKKGTPWEGGLFKLRMLFKDDYPSSPPKCKFEP
PLFHPNVYPSGTVCLSILEEDKDWRPAITIKQILLGIQELLNEPNIQDPAQAEAYTIYCQNRVEYEKRVRAQAKKFAPS

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储:+4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.