

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

UXT, 1-157aa, Human, His tag, E.coli

产品货号: TP04448

第三版

别名:Protein UXT, ART-27, STAP1

**描述:**Ubiquitously-expressed transcript, also known as UXT, belongs to the UXT family. This protein functions as a cofactor that modulates androgen receptor-dependent transcription, and also plays a critical role in tumor necrosis factor-induced apoptosis. Expression of this gene may play a role in tumorigenesis. Recombinant human UXT protein, fused to His-tag at N-terminus, was expressed in E.coli.

**配方:**Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol

**分子量:**20.4 kDa (177aa)

**序列:**

MGSSHHHHHHSSGLVPRGSHMATPPKRRAVEATGEKVLRYETFISDVLQRDLRKVLDHRDKVYEQLAKYLQLRNVIERL  
QEAKHSELYMQVDLGCNFFVDTVVPDTSRIYVALGYGFFLELTLAEALKFIDRKSSLLTELSNSLTKDSMNIKAHIIHMLLEGL  
RELQGLQNFPEKPHH

**纯度:**> 95% by HPLC

**浓度:**1 mg/ml (determined by Bradford assay)

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.