

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

YAE1D1, 1-226aa, Human, His tag, E.coli

产品货号: TP04516

第三版

别名: Yae1 domain containing 1, C7orf361, GK003

描述: Recombinant human YAE1D1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1mM DTT

分子量: 27.7 kDa (249aa), confirmed by MALDI-TOF (Molecular size on SDS-PAGE will appear higher)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMWVQAASLIQGPDKGDVFDEEADESLLAQREWQSNMQRRVKEGYRDGIDAG  
KAVTLQQGFNQGYKKGAEVILNYGRLRGTLSALLSWCHLHNNNSTLINKINLLDAVGQCEEYVLKHLKSITPPSHVVDL  
LDSIEDMDLCHVVPAAEKKIDEAKDERLCENNAEFNKNCSSKSHSGIDCSYVECCRTQEHAHSENPSPTWILEQTASLVKQL  
GLSVDVLQHLKQL

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.