

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

YIF1B, 1- 156aa, Human, His tag, E.coli

产品货号: TP04521

第三版

别名:protein YIF1B isoform 5, FinGER8

描述:YIF1b belongs to the YIF1 family. Recombinant human YIF1B protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol 0.4M Urea

分子量:19.5 kDa (179aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMHPAGLAAAAAGTPRLRKWPSKRRIPVSQPGMADPHQLFDDTSSAQSRGYGAQR
APGGLSYPAASPTPHAAFLADPVSNMAMAYGSSLAAQKELVDKNIDRFIPITKLKYYFAVDTMVYVGRKLGLLFFPYLHGD
WEVQYQQDTPVAPRFDVNAPD

纯度:> 95% by HPLC

浓度:1mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.