

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

YRS (Tyrosyl-tRNA synthetase)

产品货号: TP04523

第三版

别名:

描述: Tyrosyl-tRNA synthetase belongs to the class I tRNA synthetase family. Cytokine activities have also been observed for the human tyrosyl-tRNA synthetase, after it is split into two parts, an N-terminal fragment that harbors the catalytic site and a C-terminal fragment found only in the mammalian enzyme. The Nterminal fragment is an interleukin-8-like cytokine, whereas the released C-terminal fragment is an EMAP II-like cytokine.

配方: 1X PBS(pH7.2) 200mM NaCl 20% Glycerol

分子量: 59kDa

序列:

纯度: > 95% by HPLC

浓度: 0.50mg/ml

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.