



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

ZNF32, 1-273aa, Human, His tag, E.coli

产品货号: TP04537

第三版

别名: Zinc finger protein 32, KOX30

描述: ZNF32 belongs to the krueppel C2H2-type zinc-finger protein family and contains 7 C2H2-type zinc fingers. The protein may be involved in transcriptional regulation. Recombinant human ZNF32 protein, fused to His-tag at N-terminus, was expressed in E.coli

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.4M Urea

分子量: 33.4 kDa (296aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMFGFPTATLLDCHGRYAQNVAFFNVMTEAHHKYDHSEATGSSSWDIQNSFRREK  
LEQKSPDSKTLQEDSPGVRQRVYECQECGKSFQRQKGSLLTHERIHTGQKPFECTHCGKSFRAKGNLVTHQRIHTGEKPYQ  
CKEKGKSFQKSLAVHERLHTGQKPYECAICQRSFRNQSNLAVHRRVHSGEKPYRCDQCGKAFSQKGLIVHIRVHTGL  
KPYACTQCRKSFHTRGNCILHGKIHTGETPYLCGQCGKSFTQRGSLAVHQRSCSQRLTL

纯度: > 95% by HPLC

浓度: 1.0 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.