

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Small Ubiquitin-Related Modifier 3/SUMO3/SMT3A Protein(N-6His)

产品货号: TP04549

第三版

别名:Small Ubiquitin-Related Modifier 3; SUMO-3; SMT3 Homolog 1; SUMO-2; Ubiquitin-Like Protein SMT3B; Smt3B; SUMO3; SMT3B; SMT3H1

描述:SUMO3 belongs to the SUMO protein family and operates like ubiquitin. Ubiquitin-like protein which can be covalently attached to target lysines either as a monomer or as a lysine-linked polymer. Nevertheless unlike ubiquitin that targets proteins for degradation, SUMO3 takes part in several cellular processes, such as nuclear transport, transcription regulation, apoptosis and protein stability. SUMO3 participates in amyloid beta generation and has a key role in the onset or progression of Alzheimer's disease.

配方:PBS

分子量:13.8 kDa

序列:Met 1-Phe103

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.