

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Histone Chaperone ASF1A/ASF1A Protein(C-6His, N-T7 tag)

产品货号: TP04567

第三版

别名: Histone Chaperone ASF1A; Anti-Silencing Function Protein 1 Homolog A; hAsf1; hAsf1a; CCG1-Interacting Factor A; CIA; hCIA; ASF1A

描述: Human Histone Chaperone ASF1A (ASF1A) belongs to the H3/H4 family of histone chaperone proteins. ASF1A is ubiquitously expressed in many cells and tissues, interacting with histones H3 and H4. ASF1A cooperates with Chromatin Assembly Factor 1 to promote replication-dependent chromatin assembly and with HIRA to promote replication-independent chromatin assembly. In addition, ASF1A is necessary for the formation of senescence-associated heterochromatin foci (SAHF) and efficient senescence-associated cell cycle exit.

配方: PBS

分子量: 25.4 kDa

序列: Met 1-Met204

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.