

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Ubiquitin-Conjugating Enzyme E2 H Protein

产品货号: TP04574

第三版

别名: Ubiquitin-Conjugating Enzyme E2 H; UbcH2; Ubiquitin Carrier Protein H; Ubiquitin-Conjugating Enzyme E2-20K; Ubiquitin-Protein Ligase H; UBE2H

描述: Ubiquitin-Conjugating Enzyme E2 H (UBE2H) belongs to the E2 Ubiquitin-Conjugating Enzyme family. The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. It has been shown to conjugate ubiquitin to histone H2A in an E3 dependent manner in vitro. UBE2H is the human homolog to the yeast DNA repair gene RAD6, which is induced by DNA damaging reagents. UBE2H has been associated with cancer-induced cachexia and with the regulation of sepsis-induced muscle proteolysis.

配方: PBS

分子量: 47.0 kDa

序列: Met 1-Leu183

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.