

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Wnt Inhibitory Factor 1/WIF-1 Protein(C-6His)

产品货号: TP04580

第三版

别名:Wnt Inhibitory Factor 1; WIF-1; WIF1

描述:Wnt Inhibitory Factor 1 (WIF1) is a secreted protein, which binds WNT proteins and inhibits their activities. WNT proteins are extracellular signaling molecules involved in the control of embryonic development. WIF1 contains a WNT inhibitory factor (WIF) domain and 5 epidermal growth factor (EGF)-like domains. is found to be present in fish, amphibia and mammals. WIF1 is a recurrent target in human salivary gland oncogenesis.WIF1 may be involved in mesoderm segmentation. WIF1 is a tumor suppressor, specifically in nonfunctioning pituitary tumors.

配方:PBS

分子量:39.5 kDa

序列:Gly29-Trp379

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.