



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Sentrin-Specific Protease 8/SEN8 Protein(N-6His)

产品货号: TP04594

第三版

别名: Sentrin-Specific Protease 8; Deneddylase-1; NEDD8-Specific Protease 1; Protease Cysteine 2; Sentrin/SUMO-Specific Protease SENP8; SENP8; DEN1; NEDP1; PRSC2

描述: Sentrin-Specific Protease 8 (SEN8) mediates the reversible covalent modification of proteins by NEDD8. SEN8 catalyzes the full-length NEDD8 to generate its mature form and deconjugation of NEDD8 from targeted proteins such as CUL2 , CUL4A in vivo, or p53. but it does not show activity against ubiquitin or SUMO proteins.

配方: PBS

分子量: 26.2 kDa

序列: Met 1-Lys212

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.