

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Death Receptor 6 Protein

产品货号: TP04612

第三版

别名:Tumor Necrosis Factor Receptor Superfamily Member 21; Death Receptor 6; CD358; TNFRSF21; DR6

描述:Tumor Necrosis Factor Receptor Superfamily Member 21 (TNFRSF21) is a type I transmembrane receptor that includes four extracellular cysteine-rich motifs and a cytoplasmic death domain. DR6 is highly expressed in heart, brain, placenta, pancreas, lymph node, thymus and prostate. DR6 may activate NF-kappa-B and JNK to promote apoptosis and T-cell differentiation. In addition, DR6 binds with N-APP, which is released by the deprivation of Trophic-factor. It triggers caspase activation and degeneration of both neuronal cell bodies (via caspase-3) and axons (via caspase-6). DR6 is also expressed on the tumor cell lines and can be induced by TNF- α .

配方:PBS

分子量:34.6 kDa

序列:Gln42-Leu350

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.