

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human IL-18 Receptor Accessory Protein/IL-18RAcP/IL-1R7 Protein(C-6His)

产品货号: TP04617

第三版

**别名:** Interleukin-18 receptor accessory protein (IL18RAP); Accessory protein-like; CD218 antigen-like family member B; CDw218b; IL-1R accessory protein-like; Interleukin-1 receptor 7; Interleukin-18 receptor accessory protein-like; Interleukin-18 receptor beta  $\epsilon$  " CD218b and IL1R7.

**描述:** IL-18RAP is a single-pass type I membrane protein and contains two Ig-like C2-type domains and one TIR domain, IL18RAP can be induced by IFN-alpha and IL12 in nature killer cells and T-cells. The coexpression of IL18R1 and IL18RAP is required for the activation of NF-kappa B and JNK in response to IL-18. Defects in IL18RAP are associated with Coeliac disease.

**配方:** PBS

**分子量:** 39.4 kDa

**序列:** Phe20-Arg356

**纯度:** > 95% by HPLC

**浓度:**

**内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.