

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Complement Component C1q Receptor/C1QR1/CD93 Protein(C-6His)

产品货号: TP04644

第三版

别名: Complement Component C1q Receptor; C1q/MBL/SPA Receptor; C1qR; C1qR(p); C1qRp; CDw93; Complement Component 1 q Subcomponent Receptor 1; Matrix-Remodeling-Associated Protein 4CD93; C1QR1; MXRA4

描述: C1q R1 is also known as CD93, collectin receptor, and AA4 antigen, belongs to the Group XIV C-Type lectin family which play a role not only in cell – cell adhesion processes but also in host defence. All of them contain a C-type lectin domain, a series of epidermal growth factor like domains (EGF), a highly glycosylated mucin-like domain, a unique transmembrane domain and a short cytoplasmic tail. C1q R1 has also been identified as a hematopoietic stem cell marker.

配方: PBS

分子量: 59.2 kDa

序列: Ala24-Lys580

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.