

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Secretogranin-3 Protein

产品货号: TP04653

第三版

别名: Secretogranin-3; Secretogranin III; SgIII;

描述: Secretogranin-3 (SCG3) is a member of the chromogranin/secretogranin family of neuroendocrine secretory proteins. SCG3 is expressed in brain, heart, kidney, liver and skeletal muscle. It may serve as precursors for biologically active peptides. Some granins have been shown to function as helper proteins in sorting and proteolytic processing of prohormones; however, the function of this protein is unknown.

配方: PBS

分子量: 52.0 kDa

序列: Phe20-Leu468

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.